

47 Solvents in working environment air

Column InertCap 1
 0.53mm I.D. × 30m df=2.0µm
 Column Temp. 40°C (5min hold) → (4°C / min) → 200°C (10min hold)
 Injector Temp. 250°C
 Detector Temp. 250°C
 Injection Size 0.2µL

1 Methanol	14 i-Butanol	24 N,N-Dimethylformamide	35 Styrene
2 Ethanol	Methylcellosolve	25 n-Amyl alcohol	36 Cellosolve acetate
3 Acetone	Tetrahydrofuran	26 Toluene	37 o-Xylene
4 i-Propanol	15 1,2-Dichloroethane	27 i-Butyl acetate	1,1,2,2-Tetrachloroethane
5 Ethyl ether	16 1,1,1-Trichloroethane	28 Methyl butyl ketone	38 Butylcellosolve
6 Methyl acetate	17 i-Propyl acetate	29 n-Butyl acetate	39 n-Amyl acetate
7 Dichloromethane	n-Butanol	30 Tetrachloroethylene	40 Methylcyclohexanone
8 n-Propanol	18 Benzene	31 Chlorobenzene	41 Methylcyclohexanol
Carbon disulfide	19 Carbontetrachloride	32 Ethylbenzene	42 Phenol
9 1,2-Dichloroethylene	20 Trichloroethylene	33 m-Xylene	43 o-Dichlorobenzene
10 Methyl ethyl ketone	1,4-Dioxane	p-Xylene	44 o-Cresol
11 2-Butanol	21 Ethylcellosolve	i-Amyl acetate	45 m-Cresol
12 Ethyl acetate	22 n-Propyl acetate	34 Cyclohexanone	p-Cresol
13 Hexane	23 i-Amyl alcohol	Cyclohexanol	
Chloroform	Methyl i-butyl ketone		

