




Evolution

Workstation Software for OPTIC 3 Injector


The newly designed, highly innovative software expands the dimensions of injection technology. Based on years of experience it extends OPTIC 3 features and optimises the analytical output. Evolution Workstation offers state of the art GC injection control in a user-friendly way. Integrated Internet technology secures support and software updates at any time at any location.

Features


- Complete status information at any moment during standby state or method run
- Real-time graphical run-time parameters display
- Easy, on-click analytical method definition and development
- Automatic generation of a method optimisation sequence
- Possibility to save the run-time data for every injection
- Built-in At Column calculator for most used solvents
- Direct control of the instrument in standby state
- System and method log files
- Password protection with two access levels
- On-line help
- RS 232 or LAN interface to the instrument
- Windows 98SE/2000/NT/XP compatibility



Method development



Optimization Sequence



Captured explosion

