

Indoor air analysis by AutoCan

Technique GC
Column AQUATIC
0.25mm I.D. x 60m df=1.0µm

- | | |
|---------------------------|-----------------|
| 1 Propane | 8 Ethylacetate |
| 2 Unknown | 9 Toluene |
| 3 Butane | 10 Ethylbenzene |
| 4 Methanol / Acetaldehyde | 11 m-Xylene |
| 5 2-Methylbutane | p-Xylene |
| 6 Ethanol | 12 α-Pinene |
| 7 Acetone | 13 Limonene |

